

Automation Environment Setup

We need to perform the following process to set Automation Environment in our machine.


Installing Python:

1. We need to download python version 3.8.5 from the following link.

<https://www.python.org/downloads/>

Note: For now(**Milestone-6**) we have installed **python version 3.8.5**


2. While installing python we need to **check 'Add Python 3.8 to Path'** checkbox from Python Installation Window.
3. We need to make sure to choose custom installation, this this we will create a custom directory "C:/Python" and install python there to keep the directory path short and simple.



With this, Python 3.8 will be added to Environment variables of the computer.

Installing Source Tree and Clone Repository:


1. We need to download source tree from the following link and install that in our machine.
<https://www.sourcetreeapp.com/>
2. After installation is complete we have to clone the master branch in our machine. For doing that we need to do the following.
 - Navigate to https://dev.azure.com/planetdds/PoWProjects/_git/PoWTestCases
 - Copy the Repository URL



- Paste the Repo URL in Source Path Field of Clone Section of Source Tree.


- **Be careful not to change mouse focus outside of the input field in this stage!** Replace the 'planetdds' with the username (mail prefix) of PlanetDDS email address (ex: **username**@contract.planetdds.com) And select the location where we want to clone the repo (Preferably "D://PowTestCases").

3. Click on "Advanced Options". Then change the checkout branch to "Main". Then click the clone button.



Installing DiffMerge and Integrate with Source Tree:

1. Download DiffMerge <https://sourcegear.com/diffmerge/downloads.php> and install this to the machine
2. After completing installation open Source Tree and go to Tools->Options
3. Switch to Diff Tab and Set External Diff Tool and Merge Tool to DiffMerge and Click OK button.



Installing Python Libraries:

1. For installing python libraries we run a bat file which can be found in cloned repository. In our cloned repository **environment_setup.bat** can be found on the following directory.
'DenticonTestAutomation\DenticonTestAutomation\TestAutomation'
2. Run the **environment_setup.bat**

Local Disk (D:) > Auto > DenticonTestAutomation > DenticonTestAutomation > TestAutomation				
Name	Date modified	Type	Size	
.idea	2/26/2021 2:56 PM	File folder		
RobotFramework	2/15/2021 1:19 PM	File folder		
TestFramework	2/25/2021 7:22 PM	File folder		
_Denticon_Test_Runner	2/26/2021 2:35 PM	Executable Jar File	40 KB	
day_in_life_ci_cd_runner_dev	2/15/2021 1:19 PM	Windows Batch File	1 KB	
day_in_life_ci_cd_runner_production	2/15/2021 1:19 PM	Windows Batch File	1 KB	
day_in_life_runner_beta	2/26/2021 2:38 PM	Windows Batch File	2 KB	
day_in_life_runner_dev	2/15/2021 1:19 PM	Windows Batch File	2 KB	
day_in_life_runner_production	2/26/2021 2:38 PM	Windows Batch File	2 KB	
day_in_life_runner_sandbox	2/26/2021 2:38 PM	Windows Batch File	2 KB	
day_in_life_runner_staging	2/26/2021 2:38 PM	Windows Batch File	2 KB	
environment_setup	2/26/2021 2:43 PM	Windows Batch File	1 KB	
prerequisite_test_runner	2/15/2021 1:19 PM	Windows Batch File	2 KB	
regression_test_runner_production	2/26/2021 2:38 PM	Windows Batch File	2 KB	
regression_test_runner_staging	2/26/2021 2:38 PM	Windows Batch File	2 KB	
test_runner	2/26/2021 2:38 PM	Windows Batch File	2 KB	

Installing JDK and Running Denticon Test Runner:

1. Install the latest JDK from the link below. Please be sure to install only the stable releases.

[Download the Latest Java LTS Free](#)

2. We can now run test automation from an user friendly UI we can use 'Denticon Test Runner'. In our cloned repository can be found on the following directory.

'DenticonTestAutomation\DenticonTestAutomation\TestAutomation'


Local Disk (D:) > Auto > DenticonTestAutomation > DenticonTestAutomation > TestAutomation				
Name	Date modified	Type	Size	
.idea	2/26/2021 2:56 PM	File folder		
RobotFramework	2/15/2021 1:19 PM	File folder		
TestFramework	2/25/2021 7:22 PM	File folder		
_Denticon_Test_Runner	2/26/2021 2:35 PM	Executable Jar File	40 KB	
day_in_life_ci_cd_runner_dev	2/15/2021 1:19 PM	Windows Batch File	1 KB	
day_in_life_ci_cd_runner_production	2/15/2021 1:19 PM	Windows Batch File	1 KB	
day_in_life_runner_beta	2/26/2021 2:38 PM	Windows Batch File	2 KB	
day_in_life_runner_dev	2/15/2021 1:19 PM	Windows Batch File	2 KB	
day_in_life_runner_production	2/26/2021 2:38 PM	Windows Batch File	2 KB	
day_in_life_runner_sandbox	2/26/2021 2:38 PM	Windows Batch File	2 KB	
day_in_life_runner_staging	2/26/2021 2:38 PM	Windows Batch File	2 KB	
environment_setup	2/26/2021 2:43 PM	Windows Batch File	1 KB	
prerequisite_test_runner	2/15/2021 1:19 PM	Windows Batch File	2 KB	
regression_test_runner_production	2/26/2021 2:38 PM	Windows Batch File	2 KB	
regression_test_runner_staging	2/26/2021 2:38 PM	Windows Batch File	2 KB	
test_runner	2/26/2021 2:38 PM	Windows Batch File	2 KB	

IDE(Pycharm) Installation and Configuration:

1. Pycharm is used as our IDE for developing automation project. For downloading the Pycharm, we can navigate to the following page and **download the community version of Pycharm**.


<https://www.jetbrains.com/pycharm/download/>

Note: In PyCharm, we open and work on the 'TestAutomation' folder. Here, we can find all of our Resource, Testsuites, Test Methods and other necessary files.




1. Necessary Plugins:


- If we open a bat file, pycharm will suggest us to install plugin '**Batch Script Support**'. We have to install that plugin.
- We need to install "**CMD Support**" from the Plugins Marketplace in order run certain Batch files from the IDE itself.




- If we are suggested to install Git from PyCharm, then click on "Download and Install"
- From File->Settings->Plugin, we have to install IntelliBot



At this point, since we have some .robot files opened in pycharm, it will prompt us to install “**Hyper RobotFramework Support**”. We will **ignore** this prompt and will **NOT INSTALL** this plugin.



After completing installations of the necessary plugins in PyCharm we should navigate to **File** → **Invalidate Caches**, upon which the following pop up will appear. We will click “**Invalidate and Restart**”.



Installing Autolt for interacting with OS pop-ups:


Autolt is a freeware BASIC-like scripting language designed for automating the Windows GUI and general scripting. It uses a combination of simulated keystrokes, mouse movement and window/control manipulation in order to automate tasks in a way not possible or reliable with other languages

Please download and install Autolt from this [URL](#).


Setting Email Environment:

We should copy the following files and paste these files to the 'C:\Python\Lib\site-packages\robotframework_metrics' folder.

File#1:



File#2:



Installing Google Chrome and Setting ChromeDriver:

1. In our automation environment, latest version Google Chrome should be downloaded and installed.
2. In 'C:\Python\Lib\site-packages\chromedriver_autoinstaller\utils.py' file, we need to change two paths. The paths are mentioned in the below image. We need to change these paths to:

```
[ 'powershell', '-command', '$(Get-Item -Path  
"$env:PROGRAMFILES\Google\Chrome\Application\chrome.exe").VersionInfo.FileVersion'],
```

```

def get_chrome_version():
    """
    :return: the version of chrome installed on client
    """
    platform, _ = get_platform_architecture()
    if platform == 'linux':
        path = get_linux_executable_path()
        with subprocess.Popen([path, '--version'], stdout=subprocess.PIPE) as proc:
            version = proc.stdout.read().decode('utf-8').replace('Chromium', '').replace('Google Chrome', '').strip()
    elif platform == 'mac':
        process = subprocess.Popen(['/Applications/Google Chrome.app/Contents/MacOS/Google Chrome', '--version'], stdout=subprocess.PIPE)
        version = process.communicate()[0].decode('UTF-8').replace('Google Chrome', '').strip()
    elif platform == 'win':
        process = subprocess.Popen(
            ['powershell', '-command', '$(Get-Item -Path "$env:PROGRAMFILES\Google\Chrome\Application\chrome.exe").VersionInfo.FileVersion'],
            stdout=subprocess.PIPE, stderr=subprocess.DEVNULL, stdin=subprocess.DEVNULL
        )
        output = process.communicate()
        if output and output[0] and len(output[0]) > 0:
            version = output[0].decode('UTF-8').strip().split()[-1]
        else:
            process = subprocess.Popen(
                ['powershell', '-command', '$(Get-Item -Path "$env:PROGRAMFILES\Google\Chrome\Application\chrome.exe").VersionInfo.FileVersion'],
                stdout=subprocess.PIPE, stderr=subprocess.PIPE, stdin=subprocess.PIPE
            )
            version = process.communicate()[0].decode('UTF-8').strip()
    else:
        return
    return version

```

Installing Denticon X-Ray Bridge Integration:

1. From Denticon we need to download 'Denticon X-Ray Bridge Integration' from <https://denticon.com/Support/HelpAndSupport.aspx> and install this to our machine.


The screenshot shows the Denticon Support website. The main content area is titled 'Denticon X-ray Bridge Services'. It includes a list of links for 'Denticon X-ray Bridge integration', 'Chrome', and 'Firefox'. A red box highlights the 'Denticon X-ray Bridge integration' link. The sidebar on the left contains various support links, and the right sidebar contains links for 'Online Patient Registration Link', 'Appointment Link', and 'Patient Portal Link'.

After installing this, user can navigate to Utility->Launch Third Party-EzDenti

Configuration to run Debugger in Pycharm:

We use main.py file under TestAutomation->TestFramework directory for debugging python code. To run debugger, following configuration should be done:

1. Open project in Pycharm
2. Navigate to Run->Edit Configuration. A window will open. Click on the plus (+) icon. Select option 'Python' from the dropdown.
3. Now, add the following in the window



- Name: Main
- Script path: path to the main.py file
- Environment variables: PYTHONUNBUFFERED=1 (should be automatically configured)
- Python interpreter: Python 3.8 (select from the dropdown)
- Working directory: Path to the TestFramework folder
- Add content root to PYTHONPATH: checked (should be automatically configured)
- Add source root to PYTHONPATH: checked (should be automatically configured)

4. After the setup, click on 'Apply' and then 'OK'.